

OptiK User Manual

Overview

OptiK is a command-line tool for unsupervised k-mer size optimization in genome analysis. It computes k-mer frequency matrices across a specified range of k values, reduces dimensionality via truncated SVD, performs clustering, evaluates clustering quality using internal metrics, and identifies the optimal k-mer size based on a rank-based consensus of these metrics.

Installation

1. Clone the repository:

```
git clone https://github.com/yourusername/OptiK.git
cd OptiK
```

2. Create a virtual environment and install dependencies:

```
python3 -m venv venv
source venv/bin/activate
pip install -r requirements.txt
```

Alternatively, use Conda:

```
conda env create -f environment.yml
conda activate optik
```

Usage

Basic command:

```
python3 optik_v2.py -i <input_dir> -o <output_dir>
```

Required Arguments

- `-i, --input`: Directory containing genome FASTA files
- `-o, --output`: Directory where results will be stored

Optional Arguments

- `--k-range K_MIN K_MAX`: Range of k-mer sizes to evaluate (default: 3 8)
- `--clusterer`: Clustering algorithm (kmeans or agglomerative; default: kmeans)
- `--plot-umap`: Generate UMAP projections for visualization
- `--random-seed`: Seed for reproducibility (default: 42)
- `--svd-components`: Number of components for truncated SVD (default: auto)

- `--n-clusters` MIN MAX: Range of cluster numbers to test (default: 3 8)
- `--threads`: Number of threads for parallel processing (default: 1)

Output Files

- `metrics.csv`: Clustering scores for each k and cluster count
- `best_k.txt`: Selected optimal k-mer size
- `cluster_assignments_kX.csv`: Cluster assignments for optimal k
- `umap_kX.png`: UMAP plot if visualization is enabled
- `svd_kX.npy`, `matrix_kX.npy`: Intermediate data matrices

How OptiK Chooses the Optimal k

For each k:

1. Compute average Silhouette, Calinski-Harabasz, and Davies-Bouldin scores across cluster counts
2. Rank each k within each metric
3. Compute cumulative rank score
4. Select k with lowest total rank
5. Break ties by favoring higher Silhouette and lower Davies-Bouldin scores

Reproducibility

- Set `--random-seed` for deterministic output
- All configuration parameters and outputs are saved to the output directory

Example

```
python3 optik_v2.py \
-i ./data/HpGP/ \
-o ./results/ \
--k-range 3 8 \
--clusterer kmeans \
--plot-umap \
--random-seed 42
```

Notes

- Input FASTA files should be high-quality assembled genomes.
- For metagenomic or draft assemblies, preprocessing and normalization may be required.

- For best performance on large datasets, use lower k_max or activate parallel threads.

License

OptiK is released under the MIT License.

Citation

Gutierrez Escobar AJ. OptiK: An Entropy-Driven Framework for Optimal k-mer Size Selection in Comparative Genomics. bioRxiv (2024).